


BT-710 PULSE OXIMETER



4.3" Pulse Oximeter

SpO₂, Pulse

Touch screen

Handheld design / PC view software / Over 8 hours battery use

:: Technical Specification

Model		BT-710	
Category		Pulse Oximeter	
Display		4.3" Touch LCD	
SpO ₂	Measurement Range	0-100 %	
	Accuracy (70-100%)	Adult / Pediatric	±2 %
		Neonate	±3 %
	Accuracy (0-69%)	Unspecified	
	Perfusion Index	0.05-20 %	
	Pulse Rate Range (bpm)	25-250	
Battery	Type (capacity)	Li-ion(3000 mAh)	
	Run Time	5 hour	
	Charging Time	4 hour	
PC Software Interface		RJ45, SD card slot	
Warranty		2 year	



Hand held size



Protective Cover (option)

Time	SpO ₂	PR
08-07-2017 15:51	---	---
08-07-2017 15:50	97	86
08-07-2017 15:49	98	88
08-07-2017 15:48	---	---
08-01-2017 17:02	97	79
08-01-2017 16:01	97	78
08-01-2017 16:00	98	80
08-01-2017 15:59	97	77
08-01-2017 15:57	97	79
07-21-2017 01:04	---	---
07-21-2017 01:03	---	---
07-21-2017 01:02	---	---
07-21-2017 01:01	---	---
07-21-2017 01:00	---	---
07-21-2017 00:59	---	---
07-21-2017 00:58	---	---
07-21-2017 00:57	---	---
07-21-2017 00:56	---	---
07-21-2017 00:55	---	---
07-21-2017 00:54	---	---
07-21-2017 00:53	---	---
07-21-2017 00:52	---	---

1/31
X

Trend Tabular