



**Model No. 48 (Twin)**

**Technical Data**

Intensity	1,80,000 + 1,80,000 Lux
Dome Diameter	550 mm
Depth Illumination	100-150 mm
Colour Temp.	4500-5000K
Colour Rendering Index	95 Ra
Focussing	Auto Focus
LED Average Life	>50000 Hrs
Number's of LED's	48 + 48 LED's
Power Supply	220V/50Hz AC
Power	48 + 48 Watt




**Model No. 48 (Single)**



**Technical Data**

Intensity	1,80,000 Lux
Dome Diameter	550 mm
Depth Illumination	100-150 mm
Colour Temp.	4500-5000K
Colour Rendering Index	95 Ra
Focussing	Auto Focus
LED Average Life	>50000 Hrs
Number's of LED's	48 LED's
Power Supply	220V/50Hz AC
Power	48 Watt



**Model No. 48**







**Also available:**

- **Double Dome LED OT Light**  
(Solitaire 48+48)
- **Double Dome LED OT Light**  
with Camera and Monitor

**Technical Data:**

**Intensity:** 1,80,000 Lux + 1,80,000 Lux

**Shadowless Feature:** Yes

**Numbers of LED's:** 48 + 48

**Intensity Control:** 10~100 Digital Capacitive Feather Touch (Fully Remote Controlled)

**Fitting Options Available:** Ceiling 180 Degree, Ceiling 360 Degree, Mobile, Wall Mounted

**Spot Diameter:** 200-250 mm

**Depth of Light Penetration:** 150 mm

**Colour Temperature:** 3800 – 4800 K

**Colour Rendering Index:** 95 Ra

**Focusing:** Adjustable

**LED Average life:** >50000 Hrs.

**Diameter of Light:** 700 mm

**Power Supply:** 220V/50Hz AC

**Power Consumption:** 100 Watt

**Optional Special Features:** Battery Backup