



PG OTI-ONE

Electrical hydraulic operation table



- Stainless steel structure with slow rebound latex foam mattress
- Five section body layout for multiple surgical position request
- Reliable electro-hydraulic system with smooth movement & silent lifting
- Table top sliding for X-Ray and C-Arm 350mm
- Individual kidney bridge elevation system (optional)
- Built-in back up battery (optional)
- Zero function (optional)
- Override panel control on column (optional)



PG OTI-ONE

TECHNICAL FEATURES



Designed specially to meet the different demands of general surgical operations, latex foam washable mattress, concealed rubber lift locker, single remote controller.

Class I, type B

Length: 1970mm

Width: 500mm

Height adjustment: 690-1100mm

Reversed Trendelenburg: 30°

Trendelenburg: +30°

Tilt left/right: +20°/-20°

Horizontal sliding: 350mm

Net weight overall: 250Kg

Loading capacity: 250 kg

Leg plate up/down: +15°/-90°

Head rest up/down: +60°/-90°

Back plate up/down: +90°/-45°

Split leg plate: 90° (each)

Kidney bridge elevation: 120mm

Backup battery: for 80 movements (optional)


Power supply: 100-240V / 50-60Hz



Standard accessories: mattress pad (1 pc), cable controller (1 pc), arm support (1 pair), shoulder support (1 pair) body support (1 pair), leg support (1 pair), anesthesia frame (1 set), foot plate (1 pair), fixing clamp (8 pcs.), long fixing clamp (8 pcs.)



SOME OTHER ACCESSORIES




Arm Surgery Support



Body Belt



Orthopedic Traction



X-ray Board



External Kidney Bridge



Head Rest



Instrument Tray



Skull Clamp



PROGETTI S.r.l.
Strada del Rondello, 5
10028 Trofarello (Torino) - Italy
Ph. +39 011 644 738
Fax +39 011 645 822
info@progettimedical.com
www.progettimedical.com



ISO 13485:2016 CERTIFICATION



Quality system ISO certified

Distributed by:

Exclusive Distributor
Spark-Med Health Solution
27/7 Topkhana Road, Segunbagicha ,
Dhaka-1000, info@sparkmedhs.com
01799623947.