

# ATS

SURGICAL TABLE



INTUITIVE MOTION  
UNYIELDING STABILITY

# ATS

## Intuitive Motion

Built to the foremost standard of quality, with intuition and stability in mind. The ATS Surgical table exceeds the standards in patient positioning. The ATS surgical table features a C-arm compatible table top with up to a 1000 lb lift capacity. The ATS is equipped with a 180 degree rotating tabletop to maximize patient imaging access. The horizontal movement feature maximizes efficiency in patient imaging. The ATS can also be adjusted to a low height of 24 inches off from the floor to allow the surgeon better access while sitting.


## Safety

The ATS Surgical table is mobile with 4 anti-static castors. The floor lock system braces the table in place with unyielding stability. A battery backup system can continue to power the table for hours in case of a power failure. The hand controller is outfitted with safety locks to avoid unintentional table movement. Return to level can be done quickly and efficiently with the touch of a button.



# Built To The Foremost Standard Of Quality

WITH INTUITION AND STABILITY IN MIND.



**180° Rotation**



**Horizontal Slide**

- Silent electro-hydraulic movements
- Kidney Bridge
- Electronic control with safety lock
- Mobile anti-static castors with floor lock
- Horizontal Slide
- 180 degree rotation
- “Flex” positioning
- Easy access to electro-hydraulic system for service
- Optional battery backup
- Optional auto-reset / level function
- Optional wireless remote control



# ATS TECHNICAL SPECIFICATIONS:

Bed Size	2050x500mm
Height	740mm~1100mm
Reverse Trendelenburg'	Antersion $\geq$ 30°, Retroversion $\geq$ 30°
Title Left/Right	Left Oblique $\geq$ 15°, Right Oblique $\geq$ 15°
Chest Plate Adjustment	$\geq$ 25°/ $\geq$ 35°
Back Plate Adjustment	+75° ~ -20°
Head Plate Adjustment	+20° ~ -90°
Leg Plate Adjustment	$\geq$ 90°
Leg Plate Angle	0~ 90°
Rotation	0~ 180°
Longitudinal Shift	$\geq$ 300mm
Wired Remote Control	Standard
Four-level Position Reset Function	Optional
Emergency Power Supply	Optional



## Standard Accessories of ATS

Anesthesia Frame	1	Waist Supporters	2
Shoulder Supporter	2	Square Clamp II	2
Arm Rest	2	Square Clamp I	5
Leg-Crutches	2		

Exclusive Distributor  
 Spark-Med Health Solution  
 27/7 Topkhana Road, Segunbagicha ,  
 Dhaka-1000, info@sparkmedhs.com  
 01799623947.