

LED OPERATION THEATER LIGHT



Also available:
Double Dome LED OT Light

Remote Function

TOP SUN - 700/700

Technical Data:

Intensity: 2,00,000 Lux + 2,00,000 Lux

Shadowless Feature: Yes

Numbers of LED's: 48 + 48

Intensity Control: 10–100 Digital Capacitive Feather Touch (Fully Remote Controlled)

Fitting Options Available: Ceiling 180 Degree, Ceiling 360 Degree, Mobile, Wall Mounted

Spot Diameter: 200-250 mm

Depth of Light Penetration: 150 mm

Colour Temperature: 3800 – 4800 K

Colour Rendering Index: 95 Ra

Focusing: Adjustable

LED Average life: >50000 Hrs.

Diameter of Light: 700 mm

Power Supply: 220V/50Hz AC

Power Consumption: 100 Watt

Optional Special Features: Battery Backup


Yellow and White LED Combination

Sun Bright Medical Co., Ltd.



TRAVEL DIARY

📍 MEDISTAR INDIA



FEELING
PROUD

WhatsApp 05/05
+91 8851219546

